



HOFFRICHTER *aquapoint 2*

- Attachable warm air humidifier for the point 2 respiration therapy device
- Individually adjustable level of humidification from 0 to 5 according to patient needs



HOFFRICHTER *powerpackpoint 2*

- Autonomous power supply for the point 2 respiration therapy device
- Battery operation time of up to 8 hours




HOFFRICHTER *filtersystempoint 2*

- Air filtration for the point 2 respiration therapy device
- Effective filter performance ensured by dual layer filter system (coarse and fine filters)
- Allergy relief

Order Information

Description	Order No.
HOFFRICHTER <i>aquapoint 2</i> (humidifier for point 2)	00012949
HOFFRICHTER <i>powerpackpoint 2</i> (batteries for point 2)	00012846
HOFFRICHTER <i>filtersystempoint 2</i> (efficient filter system for point 2)	00012847



point 2
Dynamic | Compact | Quiet

point 2

Dynamic | Compact | Quiet

The smallest respiration therapy device in HOFFRICHTER portfolio.

The stable therapy pressure and extremely low noise level ensure safe and comfortable use. The external batteries *powerpackpoint2* enable independent electric supply (up to 8 hours mains-free operation), thus providing users with unique freedom.



HOFFRICHTER point 2 stands out for its high performance, a Flex-Line-Mode and an optimized one-click connection for its humidifier.

Features

- 3 separate FLEX levels for inspiration and expiration
- Data can be saved to an SD card
- Very low operating noise level
- Tidal volume display
- Automatic leakage and altitude compensation
- Optional extras:
 - *aquapoint 2* air humidifier
 - *powerpackpoint 2* battery (battery service life max. 8 hours)
 - *filtersystempoint 2*

Technical Data

Dimensions (W x D x H) / Weight	170 mm x 220 mm x 95 mm / 1.5 kg
Operating voltage	100 V – 240 V AC, 50 – 60 Hz or 24 V DC
Power consumption	< 12 W at 12 hPa
Battery operation (optional)	2 x 11.1 V
Pressure range	4 – 20 hPa, 4 – 25 hPa for point 2 BILEVEL ST25
Pressure variance	< 0.5 hPa without humidifier; < 0.9 hPa with humidifier
Noise level	< 30 dB(A) at 10 hPa according to DIN EN ISO 17510-1
Counter	Operating time, motor runtime, treatment time, filter age
Displays	Pressure, pressure bar, date, time, soft start ramp, alarm status*, tidal volume (VT), battery charging symbol**
Comfort functions	Automatic start/stop, mask test, soft start ramp, flex pressure settings*** I-FLEX/E-FLEX, mask alarm, automatic start after power failure, date and time display, alarm function
Other functions	VT display (On/Off), pressure unit (hPa, mbar, cmH ₂ O)

Mode / device	point 2 CPAP	point 2 AutoCPAP	point 2 BILEVEL S20	point 2 BILEVEL ST20	point 2 BILEVEL ST25
CPAP	•	•	•	•	•
APAP	–	•	–	–	–
BILEVEL S	–	–	•	•	•
BILEVEL ST	–	–	–	•	•
BILEVEL ST	–	–	–	•	•
Pressure range	4 – 20 hPa	4 – 20 hPa	4 – 20 hPa	4 – 20 hPa	4 – 25 hPa

* if activated, ** only in combination with powerpackpoint 2, *** only in CPAP and APAP modes

Scope of Supply

1 point, 1 SD-card, 1 power supply, 1 power supply cable, 1 therapy tube, 1 carrier bag, 1 manual, 1 short manual, 2 coarse filters (1pack = 2 pcs.)

Order Information

Description	Order No.
point 2 CPAP*	00012944
point 2 AutoCPAP*	00012945
point 2 BILEVEL S20*	00012970
point 2 BILEVEL ST20*	00012958
point 2 BILEVEL ST25*	00012973
<i>aquapoint 2</i> (humidifier for point 2)	00012949

* without humidifier – item can be supplied with humidifier upon request